

pipetman® **FIXED**

Big on Reliability, Low on Cost!



PIPETMAN Fixed is an air-displacement pipette, with fixed volume for maximum results and reliability. Fourteen robust models cover a large volume range from 2µL to 5000 µL, with legendary PIPETMAN accuracy and precision.

RIGHT PIPETTE

- GLP compliance - no risk of volume selection errors. Using a dedicated fixed-volume Gilson pipette for a specific test assures reliable and consistent results.

**RIGHT PRICE**

- An economical choice for clinical diagnostics, quality control, and any routine testing.

RIGHT VOLUME

- Can be readjusted by users in the laboratory, for dense or viscous fluids (once reset the new volume cannot be changed accidentally).

Trusted results

Best results guaranteed with a performance exceeding ISO 8655 specifications. PIPETMAN Fixed's specifications are more stringent than those required to be ISO 8655 compliant.

PIPETMAN F Model	Ordering reference	PIPETMAN DIAMOND Tips*		Volume (µL)	Maximum Permissible Errors			
					Gilson		ISO 8655	
					Systematic error (µL)	Random error (µL)	Systematic error (µL)	Random error (µL)
F2	F123770	D200	DF30	2	± 0.08	≤ 0.03	± 0.08	≤ 0.04
F5	F123771	D200	DF30	5	± 0.10	≤ 0.04	± 0.125	≤ 0.075
F10	F123772	D200	DF30	10	± 0.10	≤ 0.05	± 0.12	≤ 0.08
F20	F123604	D200	DF30	20	± 0.20	≤ 0.06	± 0.2	≤ 0.1
F25	F123775	D200	DF30	25	± 0.25	≤ 0.07	± 0.5	≤ 0.2
F50	F123778	D200	DF100	50	± 0.40	≤ 0.15	± 0.5	≤ 0.2
F100	F123784	D200	DF100	100	± 0.80	≤ 0.25	± 0.8	≤ 0.3
F200	F123605	D200 D300	DF200 DF300	200	± 1.60	≤ 0.30	± 1.6	≤ 0.6
F250	F123787	D1000 D1200	DF1000 DF1200	250	± 3.00	≤ 0.75	± 4.0	≤ 1.5
F300	F123788	D1000 D1200	DF1000 DF1200	300	± 3.50	≤ 0.75	± 4.0	≤ 1.5
F400	F123789	D1000 D1200	DF1000 DF1200	400	± 3.60	≤ 0.80	± 4.0	≤ 1.5
F500	F123790	D1000 D1200	DF1000 DF1200	500	± 4.00	≤ 1.00	± 4.0	≤ 1.5
F1000	F123606	D1000 D1200	DF1000 DF1200	1000	± 8.0	≤ 1.3	± 8.0	≤ 3.0
F5000	F123607	D5000		5000	± 30	≤ 8	± 40	≤ 15.0

* Recommended tip for establishing maximum permissible errors (ISO 8655 requirement). A large range of suitable PIPETMAN DIAMOND Tips is available. (See pages 19-23).

- PIPETMAN Fixed pipettes are available with a calibration certificate issued by an ISO 17025 accredited laboratory, on request. (See pages 38-39).



Pipette Kits




All You Need to Get Started, in a Single Box!

Our pipette kits now include the PIPETMAN Classic™ & PIPETMAN Neo® ranges for all of your pipetting needs.

ALL IN ONE SAVINGS

- Each kit contains two or three PIPETMAN Classic or PIPETMAN Neo models plus samples of certified-quality PIPETMAN DIAMOND® Tips and useful accessories, such as SINGLE™ Pipette Holders (except Large Volume Kit).
- Save money when compared to pipettes and accessories bought separately.

**PIPETMAN Kits contents**

Model	Ordering reference	PIPETMAN models	PIPETMAN DIAMOND Tips (racks of 96 tips)	Accessories
PIPETMAN Neo Starter Kit	F167600	P10N P100N P1000N	DFL10ST DF100ST DF1000ST	<ul style="list-style-type: none"> • 3 SINGLE™ pipette holders • 3 microtube openers THE JIMMY™ • 1 Gilson Guide to Pipetting • 1 Two-minute inspection poster • 3 pipette identification stickers 
PIPETMAN Neo PCR Kit	F167700	P20N P200N P1000N	D200 D1000	<ul style="list-style-type: none"> • 3 SINGLE pipette holders • 3 microtube openers THE JIMMY • 1 Gilson Guide to Pipetting • 1 Two-minute inspection poster • 3 pipette identification stickers 
PIPETMAN Classic Micro-Volume Kit	F167500	P2 P10 P100	Sterile FilterTips DFL10ST DF100ST	<ul style="list-style-type: none"> • 3 SINGLE Pipette Holders • 10 COLORIS™ identification clips • 1 Gilson Guide to Pipetting 
PIPETMAN Classic Starter Kit	F167300	P20 P200 P1000	D200 D1000	<ul style="list-style-type: none"> • 3 SINGLE Pipette Holders • 10 COLORIS identification clips • 1 Gilson Guide to Pipetting 
PIPETMAN Classic Large-Volume Kit	F167400	P5000 P10mL	5 x D5000 tips 2 x D10mL tips (no racks)	<ul style="list-style-type: none"> • 2 x 10 protection filters for P5000 and P10mL 